

Digital column scale with eye-level display

seca 777

- + Medical-grade, robust scale for years of use with high patient volume
- + Built-in measuring rod with 10–230 cm measuring range
- + USB interface for EMR integration solutions or PC connectivity
- + Large rubber coated castors for excellent movability
- + Low-profile solid metal platform with tilt-proof design for increased safety
- + Display swappable to left or right side
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children



seca 777

Digital column scale with eye-level display



Absolute sturdy construction

The platform gives patients (especially those affected by obesity) a greater sense of security when stepping on it. The measuring rod is firmly attached to the column and provides a measuring range that is suitable for pediatrics and for very tall individuals. In addition, there is enough space between the platform and column so that patients can be weighed comfortably.

Cost-effective without any compromises in quality

The seca 777 puts unstable and inaccurate weighing and measuring to an end. It boasts a tilt-proof full metal platform and a stable single piece column. The eye-level display and ultra sturdy base provide comfort and security. It is a quality every patient will notice, as soon they sets foot on the platform. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.



Scale

Capacity	250 kg / 550 lbs
Graduation	100 g / 0,2 lb
Dimensions (WxHxD)	41,0 x 135,6 x 57,6 cm / 16,1 x 53,4 x 22,7 inch
Platform (WxHxD)	35,0 x 5,9 x 41,8 cm / 13,8 x 2,3 x 16,5 inch
Measuring range	10 – 230 cm / 4 – 90 inch
Weight	12,7 kg / 28,0 lbs
Power supply	Batteries, power adapter (optional)
Data transmission	USB interface for customer specific solution
Functions	TARE, HOLD, kg/lbs switch-over lock, BMI
Electromagnetic radiation resistance	60601-1-2 compliant
Accessories	+ Handrail seca 477 + Panda seca 459 + Power adapter 401

Handrail seca 477 (optional)

Dimensions (WxHxD)	42,9 x 112,2 x 30,9 cm / 16,9 x 44,2 x 12,2 inch
Weight	4,0 kg / 8,8 lbs

Quality you can trust

As the world market leader in medical measuring and weighing, we put all our experience and knowledge in the development and design of the seca 777. This includes extensive and stringent endurance tests that simulate more than ten years of daily usage. We made sure, that even after thousands of height and weight measurements and thousands of button presses, the seca 777 will work smooth and reliable as it did on the first day.