

Hospital & Primary Care





Riester

Founded in 1948 in Germany, Riester manufactures diagnostic instruments to meet the highest demands. The diagnostic products manufactured by Riester in Germany include a wide range of sphygmomanometers, otoscopes and ophthalmoscopes, and stethoscopes - to name a few.

Designed for demanding everyday use in any health care facility - from doctors' practices to hospitals, many Riester products have been created through close cooperation between their experienced research and development department and practitioners from government research centres and university hospitals.

Utilising the latest manufacturing techniques, supported by a focused investment in product development, Riester is able to offer customers the ultimate in product value in the form of high tech functionality, ease of operation, reliability and durability at affordable costs.



Diagnostic Sets

e-scope - pocket instrument 2.0

White Light, Green technology.

LED

3.7V LED illumination and conventional alkaline batteries – an almost impossible combination up until now. State of the art, the new e-scope is the perfect blend of leading-edge LED technology, more efficient diagnosis, and an environment-friendly energy source, and almost unlimited service life..

LED technology for otoscopes and ophthalmoscope with low-reflection, bright and white light. The pure white light (approx. 5500K) ensures colour-neutral illumination in addition to improve colour differentiation, thereby aiding more efficient diagnosis in comparison with halogen or xenon bulbs. The LED illumination shows colours as they really are.



IPC - Integrated Power Converter

Riester Innovation. The tiny IPC found inside the e-scope, converts the voltage of a conventional AA size battery into the optimum voltage supply for LED, thereby allowing efficient operation of LED illumination with alkaline batteries. As a result the IPC not only allows optimum light quality, but also provides economic and ecological benefits.



	Black Code	White Code
e-scope® otoscope, direct illumination, vacuum 2.7 V, in pouch	RI.2101-200	RI.2100-200
e-scope® F.O. otoscope LED 3.7 V, in case	RI.2111-203	RI.2110-203
e-scope® ophthalmoscope vacuum 2.7 V, in pouch	RI.2121-200	RI.2120-200
e-scope® ophthalmoscope LED 3.7 V, in case	RI.2123-203	RI.2122-203
e-scope® otoscope, direct illumination/ophthalmoscope, vacuum 2.7V, in pouch	RI.2131-200	RI.2130-200
e-scope® F.O. otoscope / ophthalmoscope LED 3.7 V, in case	RI.2131-203	RI.2130-203

Diagnostic


ri-scope L Otoscope

High-end Version L3 with operation lens with 4-fold magnification and specula ejection device


- glass optics with high depth of focus and precise 3-times magnification up to the image margin
- reflection minimised, consistent illumination
- fiberoptics provide highest light efficiency



Rheotronic



Reflection-reduced, cost-effective direct illumination with 2.5-V HL halogen or 3.5-V xenon lamps.



Fiber optics with 21.5-V/3.5-V XL xenon or 3.5-V LED illumination ensure an unhindered, glare-free view.



Rheotronic® makes examinations easier than flipping a switch. Riester handles with the patented rheotronic® function switch on automatically to the touch and shut off after 120 seconds by themselves when put away. The light intensity can be adjusted with the rheotronic®.

ri-scope L Ophthalmoscope

The ri-scope L ophthalmoscope's with high-performance optics, Available with correction lenses with up to 89 diopter values. High-end version L3 includes a focussing device with zero retention for force for the quick correction of ametropias.

- Dust-proof
- aspherical condenser lens
- coaxial observational and illuminational light beams
- reflection-reduced view even with small pupils.

Ri-Scope L2: Gap, fixation star, larger circle, medium circle, small circle, semi-circle. Red-free filter, blue filter, polarisation filter.



Ri-Scope L3: Grid, gap, fixation star, larger circle, medium circle, small circle, semi-circle. Red-free filter, blue filter, polarisation filter.




Diagnostic

LED + Li-Ion
inside

LED

For uniform illumination of the examination field with pure white light (550K). The LED illumination shows colours as they really are.


Helpful for efficient diagnostics plus an almost unlimited service life.



Lithium-ion (Li-Ion) battery

Memory effect and self-discharge are no match: the light-weight lithium-ion battery used in the ri-scope L has an unsurpassed lifetime - with high loading capacity and sustained high-power performance, even after well over 1000 charge cycles.

More light-weight, more powerful. Rheotronic- electronic light intensity adjustment for LEDs and light bulbs.



ri-scope L sets

Premium hand held LED diagnostic instruments. LED illumination – bright white light with a nearly unlimited lifetime!

Available with Lithium Ion rechargeable batteries for ri-charger, or with individual Lithium batteries, with a lifetime of over 100 hours of use.

ri-scope® L otoscope L1 HL 2.5 V, C-handle for 2 alkaline-batteries C-type or ri-accu®	RI.3700
ri-scope® L F.O. otoscope L2 LED 3.5 V, C-handle for 2 lithium-batteries	RI.3703
ri-scope® L F.O. otoscope L2 LED 3.5 V, C-handle for ri-accu® L	RI.3703-550
ri-scope® L otoscope L1 HL 2.5 V, AA-handle for 2 alkaline-batteries AA-type or ri-accu®	RI.3705
ri-scope® L ophthalmoscope L2 HL 2.5 V, C-handle for 2 alkaline-Batteries C-type or ri-accu®	RI.3723
ri-scope® L ophthalmoscope L2 XL 3.5 V, AA-handle for ri-accu®L	RI.3729-550
ri-scope® L ophthalmoscope L2 LED 3.5 V, C handle for 2 x Lithium batteries	RI.3725-203
ri-scope® L F.O. oto- L3/ophthal. L2 XL 2.5 V,C-handle for 2 alkaline-batteries C-type or ri-accu®	RI.3745
ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with 2 Lithium Batteries	RI.3746-203.004
ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with ri-accu®L charger (one handle, one battery)	RI.3817-203.003
ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with ri-accu®L charger (two handles, two batteries)	RI.3817-203.011



Diagnostic Interchangeable heads

LED


ri-scope L ri-derma



Proven epiluminescence microscopy allows melanocytic and non-melanocytic structures to be distinguished and enables precise observation of changes in skin pigmentation.

Illumination of the examination field with 3.5-V LED illumination resembles natural daylight. 10-fold magnification of the structure of the skin with focussing, achromatic lens.

Two different skin-friendly and sterilisable contact plates are available - With scale from 0-10mm for exact measurement of pigmented skin lesions, and without scale.



ri-derma® allows for better assessment of the skin surface structure, color and pigmentation, making it easier to distinguish between melanocytic and non-melanocytic skin changes.

ri-derma® LED 3.5 V, C handle for 2 lithium-batteries	RI.3777
ri-derma® dermatoscope head LED 3.5 V, anti-theft	RI.10577-301
Dermatoscopy oil 4 x 12.5 ml bottles	RI.10970

ri-scope L operation lens



- Small in diameter for insertion of instruments with pin connector
- 4-fold magnification
- Optical glass

ri-scope® operational lens	RI.10561-301
----------------------------	--------------

ri-scope L F.O. tongue blade holder



- 3.5V LED, XL 3.5V or XL 2.5V xenon lamp.
- Slide attachment for fixing and detaching wooden tongue depressors
- Light casing made of ABS plastic for long service life

ri-scope® F.O. tongue blade holder LED 3.5 V, anti-theft	RI.10574-301
--	--------------

ri-scope L Nasal speculum



- Swivelling and removable lense with 2.5 fold magnification
- Impact resistant casing made of glass-fibre reinforced plastic.
- 3.5V LED, XL 2.5V/3.5 V xenon lamp.

ri-scope® F.O. nasal speculum LED 3.5V, anti-theft	RI.10575-301
--	--------------

ri-scope L F.O. light conductor



- 3.5V LED, XL 2.5V/3.5 V xenon lamp.
- Best illumination of throat and laryngeal areas
- Chrome-plated metal, durable and sturdy

ri-scope® F.O. light conductor LED 3.5V anti-theft	RI.10576-301
--	--------------

ri-scope L retinoscope



- Xenon lamp.
- slit-light with a linear light beamDust-proof, very sturdy and light casing made of impact-resistant ABS plastic.

ri-scope® retinoscope head slit XL 3.5 V, anti-theft	RI.10544-301
--	--------------

ri-scope L otoscope



- 3.5V LED, HL 2.5V halogen, or XL 3.5 V xenon lamp.
- Optical Glass, reduced reflection
- Continuously variable light adjustment with auto shutdown after 120 seconds.

ri-scope®L otoscope head LED 3.5 V, anti-theft	RI.10565-301
--	--------------

ri-scope L ophthalmoscope



- HL 2.5V halogen, or XL 3.5 V xenon lamp, L2, L3.
- High performance optics, aspherical condenser lens and reduced reflection.

ri-scope®L ophthalmoscope head XL 3.5 V, anti-theft	RI.10571-301
---	--------------

Diagnostic ri-former


LED

Modular LED diagnostic station. Ri-former LED – modular wall station for ENT instruments, ophthalmic instruments, dermatoscope, sphygmomanometer, and thermometer.

The basic module is available with one handle for one instrument and can be supplemented with up to 4 more handle modules, a sphygmomanometer and a thermometer. With its centralised power supply and LED illumination, ri-former LED is practically maintenance free.



Module Combination Sample



3 year warranty for LED version, including the LEDs. White, bright, and cool LED illumination.

The handles switch on and off automatically when removed or replaced in the handle unit. The adjustable light intensity setting is maintained until next examination. Ergonomic and non-slip.

Rheotronic: enables continuously adjustable brightness as well as manual switching on and off. The spiral cord is 3m long.

ri-former wall model - Includes (1) LED (without clock);	PMD-DIA-0863
1x RI.3652-300 Wall diagnostic ri-former® 2 handles 3.5 V / 230 V without clock	
1x RI.10565-301 ri-scope® L2 otoscope head LED 3.5 V, anti-theft	
1x Ri.10595-301 Ri-scope® L2 ophthalmoscope head LED 3.5 V, anti-theft	
ri-former wall model - Includes (1) LED (with clock);	PMD-DIA-0864
1x RI.3650-300 Wall diagnostic ri-former® 2 handles 3.5 V / 230 V with clock	
1x RI.10565-301 ri-scope® L2 otoscope head LED 3.5 V, anti-theft	
1x Ri.10595-301 Ri-scope® L2 ophthalmoscope head LED 3.5 V, anti-theft	
ri-former with mobile base model - Includes (1) LED;	PMD-DIA-1181
Extension Modules	
Extension module 1 handle ri-former (2)	RI.300
Extension module ri-thermo® N ri-former® (3)	RI.3656
Extension module ri-spec ear specula dispenser (4)	RI.3654
Extension module big ben® square, adult size velcro cuff, ri-former® (5)	RI.3655-109LF
ri-former with mobile base model - Includes (1) LED;	PMD-DIA-1181
1x RI.3652-500.002 Mobile diagnostic ri-former® 2 handles 3.5 V / 230 V without clock, RI.3654 specula holder (4) and RI.10385 roll stand.	
1x RI.10565-301 ri-scope® L2 otoscope head LED 3.5 V, anti-theft	
1x Ri.10595-301 Ri-scope® L2 ophthalmoscope head LED 3.5 V, anti-theft	



Sphygmomanometers

Big Ben - Aneroid




Exact precision that lasts a lifetime.

The Big Ben analogue sphygmo offers genuine Riester quality with lifetime calibration for everyday use. Available as desk, wall, floor and rail design with a complete range of velcro cuffs available.



- Lifetime measuring technology that withstands pressures of up to maximum 600 mmHg
- Specially hardened copper-beryllium membrane, which is resistant to aging.
- Precision mechanism mounted on a special support, in order to ensure complete stability of the measuring system and the zero-point
- No zero-point fixation required
- Generously dimensioned scale (Ø 147.2 mm or 145,5 mm) and high-contrast scale layout for optimum readability; scale up to 300 mm Hg
- Maximum error tolerance: ± 3 mmHg
- Available supplied with all sized cuffs.



big ben® round desk, adult size velcro cuff	RI.1453LF
big ben® round wall, adult size velcro cuff	RI.1459LF
big ben® round floor (mobile) Ø53 cm, adult size velcro cuff	RI.1468LF
big ben® round anaesthetic, adult size velcro cuff w/o universal clamp	RI.1488LF
Universal clamp for big ben anaesthetic RI.1488LF	RI.10384
Large basket for big ben floor RI.1468LF	RI.10393
Small basket for big ben floor RI.1468LF	RI.10392

Empire - Mercury

Mercury sphygmomanometer in ABS plastic casing.

An attractive design combined with sturdy everyday-use quality. Available as desk, wall, floor, or anaesthesia model. Reasonably priced, practical and safe. Attractive optics and tested everyday practicality to cover every situation in general practice. Its specially designed stand provides a very high level of stability.

- Available as desk, mobile floor, wall, anaesthetic.
- 99.99% pure mercury
- Precision air release valve, wear-free and with fine regulation
- Easily legible tilted scale up to 300 mmHg
- Mercury lock
- Special seal against mercury contamination
- Resistant plastic housing
- Cleaning device for glass column
- Maximum error tolerance of +/- 3 mmHg



empire® N desk, adult size velcro cuff	RI.1092LF
empire® N wall, adult size velcro cuff	RI.1102LF
empire® N anaesthetic, adult size velcro cuff	RI.1132LF
empire® N floor Ø53 cm, adult size velcro cuff	RI.1112LF

Sphygmo Accessories

Bulb black	RI.10361LF
Bulb with air release valve, Riester	SPH-1LF
Tube connectors part I short + part II long, chrome plated	RI.10372
Coiled tubing, black	SPH-7LF
Air release valve, chrome plated	RI.10363LF

Sphygmomanometers

Shock-proof

The new standard for aneroid blood pressure measurement.
Exact precision that lasts.

Drop-tested up to a falling height of 120 cm.

5 years warranty on calibration

Focus Green Silver 2008 Design Award
International Forum Design award, 2010




- Drop-tested up to a falling height of 120 cm.
- Scale: Ø 52 mm white with blue pointer or black with red fluorescent pointer
- Valve base and air duct with metal precision valve
- Two-component ABS casing
- Metal threaded tube connection at the casing top for fast cuff exchange and ergonomic operation
- Maximum error tolerance of +/- 3 mm Hg
- Extensive selection of various cuff sizes
- Incl. Vinyl zippered bag

R1 shock-proof® black, velcro cuff, adult size, 1-tube, latex free	RI.1250-107LF
R1 shock-proof® white, velcro cuff, adult size, 1-tube, latex free	RI.1251-107LF
Set I, R1 shock-proof® white with 3 velcro cuffs, adult size, 1-tube, latex free	RI.1255LF

Sphygmomanometer Accessories

Manufactured using highest quality materials and craftsmanship.



	2 tube bladder cuff (Big Ben)	2 tube one piece cuff (Big Ben)	1 tube bladder cuff (Shock Proof)	1 tube bladder cuff and bulb kit
BP cuff & bladder, obese	BPC2-6LF	RI.153LF	BPC1-6LF	BPCB-6LF
BP cuff & bladder, adult	BPC2-5LF	RI.151LF	BPC1-5LF	BPCB-5LF
BP cuff & bladder, small adult	BPC2-4LF	RI.161LF	BPC1-4LF	BPCB-4LF
BP cuff & bladder, child	BPC2-3LF	RI.155LF	BPC1-3LF	BPCB-3LF
BP cuff & bladder, infant	BPC2-2LF	RI.159LF	BPC1-2LF	BPCB-2LF
BP cuff & bladder, newborn	BPC2-1LF	RI.157LF	BPC1-1LF	BPCB-1LF
BP cuff & bladder, thigh	BPC2-7LF		BPC1-7LF	BPCB-7LF
BP cuff & bladder, XL thigh	BPC2-8LF		BPC1-8LF	BPCB-8LF
Tube connector	RI.10373			
Inflation bulb & Valve	SPH-1LF			
Coiled tubing	SPH-7LF			
male/female connector set	RI.10372			

Sphygmomanometers

Minimus II

- Compact, single-tube, graduated-scale sphygmomanometer
- Highly polished chromium-plated, rugged metal housing with metal ring,
- Latex bulb with chromium-plated air release valve
- Spoon-shaped grip of stainless steel
- Tempered copper-beryllium diaphragm for long life
- Metal connector for simple and fast changes of cuffs
- Diaphragm pressure-loading capacity up to 600 mmHg
- Maximum error tolerance +/- 3 mmHg
- Linear aluminum scale Ø 49 mm, easily legible up to 300 mmHg



Minimus II black, adult size velcro cuff.	RI.1312LF
---	---------------------------

E-mega

- Fast inflation of the cuff due to new ergonomic spoon design
- ifelong precision due to almost non-aging precision movement, with tempered copper-beryllium diaphragm for pressure-loading capacity up to 600 mm Hg
- Disinfectable and washable one-piece cuff
- Large manometer scale (Ø 63 mm) and optimized scale and pointer design for better readability of values
- High-polished manometer surface, easy to disinfect
- Maximum error tolerance +/- 3 mm Hg
- Microfilter protects air release valve and measuring system



e-mega® white, disinfectable one-piece cuff, adult size, 1-tube, latex free	RI.1370-150LF
e-mega® black, disinfectable one-piece cuff, adult size, 1-tube, latex free	RI.1375-150LF

ri-san

- Tempered copper-beryllium diaphragm for long life
- Easily legible scale Ø 64 from 0 to 300 mmHg
- Latex-free bulb for pumping up the cuff
- Tube connection at the top of the pressure gage for unhindered blood pressure measurement
- Suitable for right or left handed user
- Maximum error tolerance of +/- 3 mmHg



ri-san® slate grey, adult size velcro cuff	RI.1512LF
ri-san® blue, adult size velcro cuff	RI.1517LF
ri-san® saffron, adult size velcro cuff	RI.1522LF

Sphygmomanometers

ri-champion

Fully automatic digital blood pressure monitor for precise and fast measurement.

The simple-to-use ri-champion is clinically validated according to BHS-standard A/A, and offers irregular heartbeat recognition as well as 99 memories, pulse tone, date and time. It operates with AA-batteries or an optional mains adapter.

- Measurement of blood pressure and pulse
- Recognition of arrhythmia
- Large display
- Optical and acoustic pulse signal



ri-champion® N, Pull through strap velcro cuff adult size	RI.1725-145LF
Pull through strap velcro cuff adult	RI.145LF
Pull through strap velcro cuff adult obese	RI.146LF
Pull through strap velcro cuff child size	RI.147LF

Babyphon

Suitable for children, all in blue! The palm-style sphygmomanometer is made especially for paediatricians offices and hospitals.

- Fast cuff inflation using ergonomic spoon design
- Lifelong precision due to long-life movement
- Pressure-loading capacity up to 600mm Hg
- Linear scale up to 300 mm Hg
- maximum error tolerance of +/- 3mm Hg



babyphon® plastic/metal, 3 velcro cuffs, 1 tube	RI.1441LF
---	-----------

Diplomat

A mercury sphygmomanometer in a metal case. Extremely precise, constantly improved and with impressive features. The 99.99 % mercury purity guarantees maximum measuring accuracy.

- Maximum error tolerance of +/- 3 mmHg
- Precision air release valve
- Easily legible metal scale up to 300 mmHg, fine regulation
- Mercury lock - Special seal against mercury contamination
- Unbreakable long-life metal case
- Large compartment for cuff + bulb
- Cleaning device for glass column



diplomat-presameter® desk, light grey adult size velcro cuff	RI.1002LF
diplomat-presameter® desk, blue, adult size velcro cuff	RI.1010-108LF
diplomat-presameter® desk, green, adult size velcro cuff	RI.1011-108LF
diplomat-presameter® desk, saffron, adult size velcro cuff	RI.1012-108LF

ri-cardio

Durability, Reliability, Accuracy.

24 hour Ambulatory Blood Pressure Measurements

Clinically validated measuring technology according to BHS (British Hypertension Society), EHS (European Society of Hypertension), AAMI (Association for the Advancement of Medical Instrumentation)

- Lightweight and compact design
- Motion-tolerant technology reduces risk of erroneous readings
- Patented, elastic cuff for optimal patient comfort



Ri-Cardio ambulatory BPM with adult and large adult cuff	RI.1760
Cuff for ri-cardio, small adult	RI.170
Cuff for ri-cardio, adult	RI.171
Cuff for ri-cardio, large adult	RI.172
Cuff for ri-cardio, extra large adult	RI.173

ri-vital

ri-vital® spot-check monitor. Engineered for the mobile and fast recording of vital signs. Clinically validated for adults as well as newborns weighing as little as 1kg. Simultaneous measurement of different vital signs within seconds.

- Measures SYS, DIA, MAP, pulse, SpO2, Temp, plethysmogram
- Memory function for measurement readings Shock resistant
- NIBP with mode for adults and neonates. Washable and disinfectable cuffs with coiled tube, extendible up to 2 m.
- SpO2 sensors with a standard cable length of 1.2 m.
- Easy handling via membrane keyboard (disinfectable) and self-explanatory menu.
- Efficient and lightweight rechargeable Li-ion battery pack (approx. 250 measurements)
- 2 years warranty



ri-vital® NIBP, SpO2 sensor adults + velcro cuff adult size	RI.1950-107
ri-vital® NIBP, SpO2 adult., velcro cuff adult, ri-thermo® N	RI.1951-107

ri-fox

Fast and precise measurement of pulse and oxygen content of blood (%SpO2) with familiar Riester quality in just seconds. Suitable for hospitals, doctor's practices and use at home.

- Small, mobile device with integrated sensor
- Simple and fast measurement by attaching to the finger
- Bright LED illumination for good visibility
- With visual pulse indicator
- LED-display of %SpO2, perfusion and pulse
- Battery lifetime of approx. 30 hours (in continuous operation)
- Automatic switch-off after 8 seconds
- 30 months guarantee



ri-fox N, pulse oximeter	RI.1905
--------------------------	---------

Thermometers

ri-thermo® N

- ri-thermo® N designed for easy and efficient reading of temperature in hospitals and clinics
- Infrared tympanic thermometer for the measurement of body temperature.
- LCD display with backlighting for easy reading - even in the dark
- Displays up to 12 stored temperatures
- Acoustic signal when measurement is completed or fever is detected
- Clinically validated measuring technology
- Automatic switch-off
- Probe covers included.



ri-thermo® N multifunctional infrared thermometer with 25 disposable probe covers	RI.1805
ri-thermo® N professional thermometer with dispenser and 100 probe covers	RI.1806
Dispenser with 100 pcs. of disposable probe covers for ri-thermo® N	RI.12707
Pack with 100 pcs. disposable probe covers (4 x 25 pcs.) for plastic box	RI.12701
Pack with 1000 pcs. disposable probe covers (40 x 25 pcs.) for plastic box	RI.12702

ri-gital

Easy to use, robust and precise clinical thermometer made with Riester quality. Can be used for oral, rectal and under-arm measurement of body temperature. Digital measurement is particularly accurate thanks to the highly sensitive temperature sensor.



With 3-digit display, memory and signal tone

- Unbreakable and water-tight
- Measurement range: 32 – 43.9°C
- Accuracy: +/- 0.1°
- Automatic switch-off after 10 min.

ri-gital®, digital thermometer incl. button cell battery	RI.1850
Pack with 40 pcs. disposable probe covers ri-gital®	RI.12703
Pack with 1000 pcs. disposable probe covers ri-gital®	RI.12706

Percussion Hammers

Reliability and stability distinguish all of these classical percussion hammers for the exact testing of reflexes. From the simplest model to the top quality buck hammer.

taylor stainless steel	RI.5030
buck	RI.5035
dejerine with needle	RI.5038
dejerine without needle	RI.5040
troemner 180 g	RI.5045
troemner 230 g	RI.5046
babinsky	RI.5052
berliner	RI.5055
queens	RI.5058



Stethoscopes Cardiophon 2.0



Specially developed for Cardiology.

New Technology for premium sound pickup and transference, suitable for the precise auscultation of adults and children. High precision double chestpiece with innovative floating membrane.

Patent Pending



Extra soft, replaceable, pivotable ear tips for ideal sealing and greater comfort

High precision double chestpiece with innovative bilateral floating membrane




- Innovative acoustic system with extraordinary auscultation for all frequency ranges
- High-precision double chest-piece made of stainless steel.
- Innovative bilateral floating membrane (patent pending) for precise acoustics, Ø 44 mm and 32 mm, with novel non-chill rims for improved adaptation on the skin.
- Anatomically shaped binaural with internal multiple spring for adjusting the contact pressure.
- Extra-soft, replaceable and pivotable ear-tips for ideal sealing at the ear canal and enhanced wear comfort.
- The cardiophon® is supplied with two pairs of replacement ear-tips, replacement membranes, and a name plate.



Stethoscope cardiophon 2.0, black, stainless steel	RI.4240-01
Stethoscope cardiophon 2.0, blue, stainless steel	RI.4240-03

Stethoscopes

Riester delivers brand new technology with the bell on the Cardiophon 2.0 and Duplex 2.0. This latest innovation is the bilateral floating membrane (patent pending) improving the precision of acoustics and auscultation - resulting in a far superior sound quality.



Duplex 2.0

- Precision double chest-piece with innovative membrane support.
- Anatomically shaped binaural with internal multiple spring for adjusting the contact pressure.
- Extra-soft, replaceable and pivotable ear-tips for perfect sealing at the ear canal and enhanced wear comfort
- Available in five colours.


Stethoscope duplex®2.0 black, aluminium	RI.4200-01
Stethoscope duplex®2.0 white, aluminium	RI.4200-02
Stethoscope duplex®2.0 blue, aluminium	RI.4200-03
Stethoscope duplex®2.0 red, aluminium	RI.4200-04
Stethoscope duplex®2.0 green, aluminium	RI.4200-05
Stethoscope duplex®2.0 black edition, aluminium	RI.4201-01

Stethoscope duplex®2.0 baby, white, stainless steel	RI.4220-02
Stethoscope duplex®2.0 baby, red, stainless steel	RI.4220-04
Stethoscope duplex®2.0, baby, green, stainless steel	RI.4220-05
Stethoscope duplex®2.0 neonatal, white, stainless steel	RI.4230-02
Stethoscope duplex®2.0 neonatal, red, stainless steel	RI.4230-04
Stethoscope duplex®2.0, neonatal, green, stainless steel	RI.4230-05

Duplex

- Double-headed chest-piece made of galvanized lightweight aluminium.
- Special membrane for perfect auscultation, Ø 48 mm, Bell, Ø 36 mm, Overall length: 77 cm.
- duplex® is supplied with a pair of replacement ear-tips and a replacement membrane.
- Also available in baby, neonatal, and teaching.

Stethoscope duplex®, black, aluminium	RI.4001-01LF
Stethoscope duplex®, slate grey, aluminium	RI.4001-02LF
Stethoscope duplex®, blue, aluminium	RI.4031-03LF
Stethoscope duplex®, saffron, aluminium	RI.4035-04LF
Stethoscope duplex®, green, aluminium	RI.4033-05LF
Stethoscope duplex®, baby, blue	RI.4041
Stethoscope duplex®, neonatal, blue	RI.4051
Stethoscope duplex®, teaching, slate grey	RI.4002-02



Anaestophon

Extra-ordinarily flat aluminium chest-piece without edges, easy to push under the cuff when measuring blood pressure.

- Special membrane: Ø 28 mm, Overall length: 77 cm.
- Matt chromium-plated binaural, which can be adjusted by turning.

Stethoscope anestophon®, black, aluminium	RI.4177-01LF
Stethoscope anestophon®, slate-grey, aluminium	RI.4177-02LF
Stethoscope anestophon®, blue, aluminium	RI.4177-03LF
Stethoscope anestophon®, saffron, aluminium	RI.4177-04LF
Stethoscope anestophon®, green, aluminium	RI.4177-05LF



Lighting

ri-focus LED Headlight

Wireless high performance LED headlight, powered by Li-batteries or rechargeable batteries. Green Efficiency.

- One hand operation
- Wireless and balanced headband with built-in battery compartment
- More efficient diagnosis due to white and high performance LED (140 lumen).
- LED service life up to 50,000 hours
- Significantly less energy consumption and heat generation.

ri-focus® LED, 6 V, incl. lithium batteries	RI.6090
ri-focus® LED, 6 V, incl. rechargeable NiMH batteries and plug charger 230 V	RI.6091



ri-magic LED

High performance medical LED light for diagnosis and procedures. Available as mobile, wall and examination couch or desk design.

Green Efficiency.

- High performance LED, 6.500 Kelvin.
- Lifetime of the lamp up to 50,000 hours.
- Significantly less heating and energy loss.
- Focusing device at the front of the light head.
- Highly flexible, long special arm can be bent in all directions.
- Rheostat handle for infinitely variable regulation of light intensity.

ri-magic® HP LED desk, AC adaptor 230V/120V	RI.6230.003
ri-magic® HP LED wall, AC adaptor 230V/120V	RI.6240.003
ri-magic® HP LED Floor AC adaptor 230V/120V	RI.6252.002



e-xam

Diagnostic penlights.

Photobiologically validated high-performance pen-lights for medical diagnosis of throat and pupil. Available with LED or Xenon illumination.



e-xam black, XL 2.5V, tongue blade holder	RI.5130-01
e-xam white, XL 2.5V, tongue blade holder	RI.5130-02
e-xam black, LED 2.5V, tongue blade holder	RI.5131-01
e-xam black, LED 2.5V, tongue blade holder	RI.5131-02
2.5 V LED for e-xam (pack of 1)	RI.12320
XL 2.5V Xenon bulb for e-xam (pack of 6)	RI.11305

Fortelux-N

Diagnostic penlights.

- 2.2 V vacuum lamp
- AAA batteries included
- Sterilisable aluminium casing.
- Available in assorted colours.



fortelux® N vacuum 2.2 V, black	RI.5075-526
fortelux® N vacuum 2.2 V, grey	RI.5077-526
fortelux® N vacuum 2.2 V, silver	RI.5074-526

fortelux® N vacuum 2.2 V, blue	RI.5071-526
fortelux® N vacuum 2.2 V, saffron	RI.5073-526
fortelux® N vacuum 2.2 V, green	RI.5072-526

Gynaecological Devices

Speculight

The speculight provides for the right illumination – and is the optimal addition to all vaginal specula.

- Battery handle C, chrome-plated
- For 2 batteries type C
- 85 cm cable with clip and 2.5 V standard illumination (vacuum)
- Light intensity can be regulated at will using the rheostat
- Optimal addition for any vaginal specula.



Speculight	RI.5240
------------	---------

Cusco

- External matt finish, matt inner finish
- No reflections thanks to matt inner surface
- Made of stainless steel
- Can be sterilized up to 134°C in an autoclave



cusco, for virgins	RI.5230
cusco Fig. 1	RI.5231
cusco Fig. 2	RI.5232
cusco Fig. 3	RI.5233

Graeve

- External polished finish, matt inner finish
- No reflections thanks to matt inner surface
- Made of stainless steel
- Can be sterilized up to 134°C in an autoclave



graeve Fig. 1	RI.5201
graeve Fig. 2	RI.5202
graeve Fig. 3	RI.5203

Tuning Forks

Riester offers a large range of the most high quality tuning forks and tuning fork sets for ear specialists and neurologists. Selected materials and first class finish ensure maximum functionality and reliability, which will last for a very long time.

- Highest quality materials: nickel steel or aluminium.
- Large selection of different frequencies: from 64 Hz to 4096 Hz.
- Tuning forks with and without clamps.
- Tuning forks with and without base for different requirements in everyday hospital and office practice

Tuning fork set I	RI.5140
Tuning fork set II	RI.5141
Tuning fork set III	RI.5142



Laryngoscopes

ri-dispo

- Prevention of cross contamination
- Made of high quality Makrolon® material
- High illumination, especially in combination with LED
- Pack of 20
- No sterilisation costs
- Compatible with F.O. handles (green standard)
- Metal ball lock mechanism

ri-dispo MAC 1	RI.12290
ri-dispo MAC 2	RI.12291
ri-dispo MAC 3	RI.12292
ri-dispo MAC 4	RI.12293

ri-dispo MIL 0	RI.12294
ri-dispo MIL 1	RI.12295
ri-dispo MIL 2	RI.12296
ri-dispo MIL 3	RI.12297
ri-dispo MIL 4	RI.12298



ri-modul

Macintosh and Miller Fiber Optic blade with removable light conductor, simple to remove and refit without tools.

- Unusually high light intensity with XL + LED illumination and excellent light transmission through strong fiber bundles.
- Optimum blade design for easy intubation and improved visibility.
- Miller blade size 00 for neonates.

ri-modul MAC 0	RI.12260
ri-modul MAC 1	RI.12261
ri-modul MAC 2	RI.12262
ri-modul MAC 3	RI.12263
ri-modul MAC 4	RI.12264
ri-modul MAC 5	RI.12265

ri-modul MIL 00	RI.12274
ri-modul MIL 0	RI.12275
ri-modul MIL 1	RI.12276
ri-modul MIL 2	RI.12277
ri-modul MIL 3	RI.12278
ri-modul MIL 4	RI.12279



ri-integral

- With reliable 2.7 V vacuum light source for good illumination of the examination area.
- Quick-to-remove lamp at the front of the blade.
- Can be autoclaved at 134° C for 5 minutes.*
- Stainless steel blades with easy-to-clean blade surface.
- Innovative blade design for simple intubation and good visibility.

ri-integral MAC 0	RI.12230
ri-integral MAC 1	RI.12231
ri-integral MAC 2	RI.12232
ri-integral MAC 3	RI.12233
ri-integral MAC 4	RI.12234
ri-integral MAC 5	RI.12235

ri-integral MIL 00	RI.12240
ri-integral MIL 0	RI.12241
ri-integral MIL 1	RI.12242
ri-integral MIL 2	RI.12243
ri-integral MIL 3	RI.12244
ri-integral MIL 4	RI.12245



Handle, LED, C, 2.5V

A sturdy handle is the basis of every reliable examination.

- Solid handles made of chromium plated metal.
- Fluted handle surface for hygienic usage and a secure grip in the hand.
- Practically wear-free, standard holding fixture for the blade.
- Solid metal contact for reliable activation of illumination.

Battery handle, C type, LED, 2.5V	RI.12240
-----------------------------------	----------



Page	Description	Part No.
Diagnostic Sets, e-scope - pocket instrument 2.0		
3	e-scope® otoscope, white, direct illumination, vacuum 2.7 V, in pouch	RI.2100-200
3	e-scope® F.O. otoscope, white, LED 3.7 V, in case	RI.2110-203
3	e-scope® ophthalmoscope, white, vacuum 2.7 V, in pouch	RI.2120-200
3	e-scope® ophthalmoscope, white, LED 3.7 V, in case	RI.2122-203
3	e-scope® otoscope, white, direct illumination/ophthalmoscope, white, vacuum 2.7V, in pouch	RI.2130-200
3	e-scope® F.O. otoscope, white / ophthalmoscope, white LED 3.7 V, in case	RI.2130-203
3	e-scope® otoscope, black, direct illumination, vacuum 2.7 V, in pouch	RI.2101-200
3	e-scope® F.O. otoscope, black, LED 3.7 V, in case	RI.2111-203
3	e-scope® ophthalmoscope, black, vacuum 2.7 V, in pouch	RI.2121-200
3	e-scope® ophthalmoscope, black, LED 3.7 V, in case	RI.2123-203
3	e-scope® otoscope, black, direct illumination/ophthalmoscope, black, vacuum 2.7V, in pouch	RI.2131-200
3	e-scope® F.O. otoscope, black / ophthalmoscope, black, LED 3.7 V, in case	RI.2131-203
Diagnostic Sets, ri-scope		
5	ri-scope® L otoscope L1 HL 2.5 V, C-handle for 2 alkaline-batteries C-type or ri-accu®	RI.3700
5	ri-scope® L F.O. otoscope L2 LED 3.5 V, C-handle for 2 lithium-batteries	RI.3703
5	ri-scope® L F.O. otoscope L2 LED 3.5 V, C-handle for ri-accu® L	RI.3703-550
5	ri-scope® L otoscope L1 HL 2.5 V, AA-handle for 2 alkaline-batteries AA-type or ri-accu®	RI.3705
5	ri-scope® L ophthalmoscope L2 HL 2.5 V, C-handle for 2 alkaline-Batteries C-type or ri-accu®	RI.3723
5	ri-scope® L ophthalmoscope L2 XL 3.5 V, AA-handle for ri-accu®L	RI.3729-550
5	ri-scope® L ophthalmoscope L2 LED 3.5 V, C handle for 2 x Lithium batteries	RI.3725-203
5	ri-scope® L F.O. oto- L3/ophthal. L2 XL 2.5 V,C-handle for 2 alkaline-batteries C-type or ri-accu®	RI.3745
5	ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with 2 Lithium Batteries	RI.3746-203.004
5	ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with ri-accu®L charger (one handle, one battery)	RI.3817-203.003
5	ri-scope® L FO Oto L2/ Ophthalmoscope L2 LED 3.5 V, C handle with ri-accu®L charger (two handles, two batteries)	RI.3817-203.011
6	ri-derma® LED 3.5 V, C handle for 2 lithium-batteries	RI.3777
6	ri-derma® dermatoscope head LED 3.5 V, anti-theft	RI.10577-301
6	Dermatoscopy oil 4 x 12.5 ml bottles	RI.10970
6	ri-scope® operational lens	RI.10561-301
6	ri-scope® F.O. tongue blade holder LED 3.5 V, anti-theft	RI.10574-301
6	ri-scope® F.O. nasal speculum LED 3.5V, anti-theft	RI.10575-301
6	ri-scope® F.O. light conductor LED 3.5V anti-theft	RI.10576-301
6	ri-scope® retinoscope head slit XL 3.5 V, anti-theft	RI.10544-301
6	ri-scope®L otoscope head LED 3.5 V, anti-theft	RI.10565-301
6	ri-scope®L ophthalmoscope head XL 3.5 V, anti-theft	RI.10571-301
Diagnostic Stations, ri-former wall/mobile		
7	ri-former wall model - Includes (1) LED (without clock);	PMD-DIA-0863
7	ri-former wall model - Includes (1) LED (with clock);	PMD-DIA-0864
7	Extension Modules	
7	Extension module 1 handle ri-former (2)	RI.300
7	Extension module ri-thermo® N ri-former® (3)	RI.3656
7	Extension module ri-spec ear specula dispenser (4)	RI.3654
7	Extension module big ben® square, adult size velcro cuff, ri-former® (5)	RI.3655-109LF
7	ri-former with mobile base model - Includes (1) LED;	PMD-DIA-1181

Page	Description	Part No.
Sphygmomanometers		
8	big ben® round desk, adult size velcro cuff	RI.1453LF
8	big ben® round wall, adult size velcro cuff	RI.1459LF
8	big ben® round floor (mobile) Ø53 cm, adult size velcro cuff	RI.1468LF
8	big ben® round anaesthetic, adult size velcro cuff w/o universal clamp	RI.1488LF
8	Universal clamp for big ben anaesthetic RI.1488LF	RI.10384
8	Large basket for big ben floor RI.1468LF	RI.10393
8	Small basket for big ben floor RI.1468LF	RI.10392
8	empire® N desk, adult size velcro cuff	RI.1092LF
8	empire® N wall, adult size velcro cuff	RI.1102LF
8	empire® N anaesthetic, adult size velcro cuff	RI.1132LF
8	empire® N floor Ø53 cm, adult size velcro cuff	RI.1112LF
8	Bulb black	RI.10361LF
8	Bulb with air release valve, Riester	SPH-1LF
8	Tube connectors part I short + part II long, chrome plated	RI.10372
8	Coiled tubing, black	SPH-7LF
8	Air release valve, chrome plated	RI.10363LF

Sphygmomanometers - handheld

9	R1 shock-proof® black, velcro cuff, adult size, 1-tube, latex free	RI.1250-107LF
9	R1 shock-proof® white, velcro cuff, adult size, 1-tube, latex free	RI.1251-107LF
9	Set I, R1 shock-proof® white with 3 velcro cuffs, adult size, 1-tube, latex free	RI.1255LF
9	BP cuff & bladder, obese - 2 tube, 2 piece	BPC2-6LF
9	BP cuff & bladder, adult -2 tube, 2 piece	BPC2-5LF
9	BP cuff & bladder, small adult -2 tube, 2 piece	BPC2-4LF
9	BP cuff & bladder, child -2 tube, 2 piece	BPC2-3LF
9	BP cuff & bladder, infant - 2 tube, 2 piece	BPC2-2LF
9	BP cuff & bladder, newborn -2 tube, 2 piece	BPC2-1LF
9	BP cuff & bladder, thigh -2 tube, 2 piece	BPC2-7LF
9	BP cuff & bladder, XL thigh -2 tube, 2 piece	BPC2-8LF
9	BP cuff & bladder, obese - 2 tube, 1 piece	RI.153LF
9	BP cuff & bladder, adult -2 tube, 1 piece	RI.151LF
9	BP cuff & bladder, small adult -2 tube, 1 piece	RI.161LF
9	BP cuff & bladder, child -2 tube, 1 piece	RI.155LF
9	BP cuff & bladder, infant - 2 tube, 1 piece	RI.159LF
9	BP cuff & bladder, newborn -2 tube, 1 piece	RI.157LF
9	BP cuff & bladder, obese - 1 tube, 2 piece	BPC1-6LF
9	BP cuff & bladder, adult -1 tube, 2 piece	BPC1-5LF
9	BP cuff & bladder, small adult -1 tube, 2 piece	BPC1-4LF
9	BP cuff & bladder, child -1 tube, 2 piece	BPC1-3LF
9	BP cuff & bladder, infant - 1 tube, 2 piece	BPC1-2LF
9	BP cuff & bladder, newborn -1 tube, 2 piece	BPC1-1LF
9	BP cuff & bladder, thigh -1 tube, 2 piece	BPC1-7LF
9	BP cuff & bladder, XL thigh -1 tube, 2 piece	BPC1-8LF

Page	Description	Part No.
9	BP cuff & bladder, obese - 1 tube, bladder, cuff, bulb	BPCB-6LF
9	BP cuff & bladder, adult -1 tube, bladder, cuff, bulb	BPCB-5LF
9	BP cuff & bladder, small adult -1 tube, bladder, cuff, bulb	BPCB-4LF
9	BP cuff & bladder, child -1 tube, bladder, cuff, bulb	BPCB-3LF
9	BP cuff & bladder, infant - 1 tube, bladder, cuff, bulb	BPCB-2LF
9	BP cuff & bladder, newborn -1 tube, bladder, cuff, bulb	BPCB-1LF
9	BP cuff & bladder, thigh -1 tube, bladder, cuff, bulb	BPCB-7LF
9	BP cuff & bladder, XL thigh -1 tube, bladder, cuff, bulb	BPCB-8LF
10	Minimus II black, adult size velcro cuff	RI.1312LF
10	e-mega® white, disinfectable one-piece cuff, adult size, 1-tube, latex free	RI.1370-150LF
10	e-mega® black, disinfectable one-piece cuff, adult size, 1-tube, latex free	RI.1375-150LF
10	ri-san® slate grey, adult size velcro cuff	RI.1512LF
10	ri-san® blue, adult size velcro cuff	RI.1517LF
10	ri-san® saffron, adult size velcro cuff	RI.1522LF
11	ri-champion® N, Pull through strap velcro cuff adult size	RI.1725-145LF
11	Pull through strap velcro cuff adult	RI.145LF
11	Pull through strap velcro cuff adult obese	RI.146LF
11	Pull through strap velcro cuff child size	RI.147LF
11	babyphon® plastic/metal, 3 velcro cuffs, 1 tube	RI.1441LF
11	diplomat-presameter® desk, light grey adult size velcro cuff	RI.1002LF
11	diplomat-presameter® desk, blue, adult size velcro cuff	RI.1010-108LF
11	diplomat-presameter® desk, green, adult size velcro cuff	RI.1011-108LF
11	diplomat-presameter® desk, saffron, adult size velcro cuff	RI.1012-108LF

Ambulatory Blood Pressure Monitor

12	Ri-Cardio ambulatory BPM with adult and large adult cuff	RI.1760
12	Cuff for ri-cardio, small adult	RI.170
12	Cuff for ri-cardio, adult	RI.171
12	Cuff for ri-cardio, large adult	RI.172
12	Cuff for ri-cardio, extra large adult	RI.173
12	ri-vital® NIBP, SpO2 sensor adults + velcro cuff adult size	RI.1950-107
12	ri-vital® NIBP, SpO2 adult., velcro cuff adult, ri-thermo® N	RI.1951-107
12	ri-fox N, pulse oximeter	RI.1905

Vital Signs

12	ri-vital® NIBP, SpO2 sensor adults + velcro cuff adult size	RI.1950-107
12	ri-vital® NIBP, SpO2 adult., velcro cuff adult, ri-thermo® N	RI.1951-107

Pulse Oximeter

12	ri-fox N, pulse oximeter	RI.1905
----	--------------------------	---------

Thermometers

13	ri-thermo® N multifunctional infrared thermometer with 25 disposable probe covers	RI.1805
13	ri-thermo® N professional thermometer with dispenser and 100 probe covers	RI.1806
13	Dispenser with 100 pcs. of disposable probe covers for ri-thermo® N	RI.12707
13	Pack with 100 pcs. disposable probe covers (4 x 25 pcs.) for plastic box	RI.12701

Page	Description	Part No.
Thermometers (continued)		
13	ri-gital®, digital thermometer incl. button cell battery	RI.1850
13	Pack with 40 pcs. disposable probe covers ri-gital®	RI.12703
13	Pack with 1000 pcs. disposable probe covers ri-gital®	RI.12706
Percussion Hammers		
13	taylor stainless steel percussion hammer	RI.5030
13	buck percussion hammer	RI.5035
13	dejerine with needle percussion hammer	RI.5038
13	dejerine without needle percussion hammer	RI.5040
13	troemner 180 g percussion hammer	RI.5045
13	troemner 230 g percussion hammer	RI.5046
13	babinsky percussion hammer	RI.5052
13	berliner percussion hammer	RI.5055
13	queens percussion hammer	RI.5058
Stethoscopes		
14	Stethoscope cardiophon 2.0, black, stainless steel	RI.4240.01
14	Stethoscope cardiophon 2.0, blue, stainless steel	RI.4240-03
15	Stethoscope duplex®2.0 black, aluminium	RI.4200-01
15	Stethoscope duplex®2.0 white, aluminium	RI.4200-02
15	Stethoscope duplex®2.0 blue, aluminium	RI.4200-03
15	Stethoscope duplex®2.0 red, aluminium	RI.4200-04
15	Stethoscope duplex®2.0 green, aluminium	RI.4200-05
15	Stethoscope duplex®2.0 black edition, aluminium	RI.4201-01
15	Stethoscope duplex®2.0 baby, white, stainless steel	RI.4220-02
15	Stethoscope duplex®2.0 baby, red, stainless steel	RI.4220-04
15	Stethoscope duplex®2.0, baby, green, stainless steel	RI.4220-05
15	Stethoscope duplex®2.0 neonatal, white, stainless steel	RI.4230-02
15	Stethoscope duplex®2.0 neonatal, red, stainless steel	RI.4230-04
15	Stethoscope duplex®2.0, neonatal, green, stainless steel	RI.4230-05
15	Stethoscope duplex®, black, aluminium	RI.4001-01LF
15	Stethoscope duplex®, slate grey, aluminium	RI.4001-02LF
15	Stethoscope duplex®, blue, aluminium	RI.4031-03LF
15	Stethoscope duplex®, saffron, aluminium	RI.4035-04LF
15	Stethoscope duplex®, green, aluminium	RI.4033-05LF
15	Stethoscope duplex®, baby, blue	RI.4041
15	Stethoscope duplex®, neonatal, blue	RI.4051
15	Stethoscope duplex®, teaching, slate grey	RI.4002-02
15	Stethoscope anestophon®, black, aluminium	RI.4177-01LF
15	Stethoscope anestophon®, slate-grey, aluminium	RI.4177-02LF
15	Stethoscope anestophon®, blue, aluminium	RI.4177-03LF
15	Stethoscope anestophon®, saffron, aluminium	RI.4177-04LF
15	Stethoscope anestophon®, green, aluminium	RI.4177-05LF

Page	Description	Part No.
Medical Lighting		
16	ri-focus® LED, 6 V, incl. lithium batteries	RI.6090
16	ri-focus® LED, 6 V, incl. rechargeable NiMH batteries and plug charger 230 V	RI.6091
16	ri-magic® HP LED desk, AC adaptor 230V/120V	RI.6230.003
16	ri-magic® HP LED wall, AC adaptor 230V/120V	RI.6240.003
16	ri-magic® HP LED Floor AC adaptor 230V/120V	RI.6252.002
16	e-xam black, XL 2.5V, tongue blade holder	RI.5130-01
16	e-xam white, XL 2.5V, tongue blade holder	RI.5130-02
16	e-xam black, LED 2.5V, tongue blade holder	RI.5131-01
16	e-xam black, LED 2.5V, tongue blade holder	RI.5131-02
16	2.5 V LED for e-xam (pack of 1)	RI.12320
16	XL 2.5V Xenon bulb for e-xam (pack of 6)	RI.11305
16	fortelux® N vacuum 2.2 V, black	RI.5075-526
16	fortelux® N vacuum 2.2 V, grey	RI.5077-526
16	fortelux® N vacuum 2.2 V, silver	RI.5074-526
16	fortelux® N vacuum 2.2 V, blue	RI.5071-526
16	fortelux® N vacuum 2.2 V, saffron	RI.5073-526
16	fortelux® N vacuum 2.2 V, green	RI.5072-526
Virginal Speculum		
17	Speculight	RI.5240
17	cusco speculum, for virgins	RI.5230
17	cusco speculum Fig. 1	RI.5231
17	cusco speculum Fig. 2	RI.5232
17	cusco speculum Fig. 3	RI.5233
17	graeve speculum Fig. 1	RI.5201
17	graeve speculum Fig. 2	RI.5202
17	graeve speculum Fig. 3	RI.5203
Tuning Forks		
17	Tuning fork set I	RI.5140
17	Tuning fork set II	RI.5141
17	Tuning fork set III	RI.5142
Laryngoscopes		
18	ri-dispo MAC 1	RI.12290
18	ri-dispo MAC 2	RI.12291
18	ri-dispo MAC 3	RI.12292
18	ri-dispo MAC 4	RI.12293
18	ri-dispo MIL 0	RI.12294
18	ri-dispo MIL 1	RI.12295
18	ri-dispo MIL 2	RI.12296
18	ri-dispo MIL 3	RI.12297
18	ri-dispo MIL 4	RI.12298

Page	Description	Part No.
Laryngoscopes (continued)		
18	ri-modul MAC 0	RI.12260
18	ri-modul MAC 1	RI.12261
18	ri-modul MAC 2	RI.12262
18	ri-modul MAC 3	RI.12263
18	ri-modul MAC 4	RI.12264
18	ri-modul MAC 5	RI.12265
18	ri-modul MIL 00	RI.12274
18	ri-modul MIL 0	RI.12275
18	ri-modul MIL 1	RI.12276
18	ri-modul MIL 2	RI.12277
18	ri-modul MIL 3	RI.12278
18	ri-modul MIL 4	RI.12279
18	ri-integral MAC 0	RI.12230
18	ri-integral MAC 1	RI.12231
18	ri-integral MAC 2	RI.12232
18	ri-integral MAC 3	RI.12233
18	ri-integral MAC 4	RI.12234
18	ri-integral MAC 5	RI.12235
18	ri-integral MIL 00	RI.12240
18	ri-integral MIL 0	RI.12241
18	ri-integral MIL 1	RI.12242
18	ri-integral MIL 2	RI.12243
18	ri-integral MIL 3	RI.12244
18	ri-integral MIL 4	RI.12245



AU:



1800 67 5432



sales@ecomед.com.au



www.ecomed.com.au

Unit 4, 30 Foundry Road
Seven Hills NSW 2147

NZ:



0800 300 421



sales@ecomед.co.nz



www.ecomed.co.nz