

LUXIFLEX LED PLUS

Luxiflex LED Plus is a cold light luminaire, versatile with a compact design. Luxiflex is designed for examination procedures that require great performance.

Flexibility

Flexible arm allows to adjust the height to assure the proper working distance.

Compact design

Luxiflex offers easy use, thanks to its compact LED head with a diameter of 8 cm and its switch located at the base of the lamp. These features make this model an ideal luminaire for the most demanding examinations procedures.

Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

High performance LED

Luxiflex LED Plus 35.000 lux at 50 centimeters.

LUXIFLEX LED PLUS

Technical Details	(Cod. I0104NL)
Light source	LED 10W 14°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	35.000 lux
Colour temperature	4.500 °K
Colour rendering	97%
Total length	85 cm
Colour	White





Wall support
Cod. B0101800



AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026



Extension wall arm
Cod. B0101S0



Trolley stand
4 kg Cod. 09655
8,8 kg Cod. 09651