

# seca 952

## Chair scale for weighing while seated

new



- The simple solution for disabled people.
- Fold-up armrests and footrests make sitting down easy.
- Brakes on wheels for secure stance.
- Intuitive operating elements simplify use.
- All essential functions for maximum usage.

**seca 952:**Relief for patient and  
medical personnel.

**S**tanding on a scale is a problem for physically disabled or older people. Chair scales are of great assistance to patients and medical staff in nursing homes, hospitals, geriatric and orthopedic departments. The new seca 952 offers a sound, user-oriented solution with all the functions that a good chair scale should offer.

## Everything that counts.

With a capacity of 200 kg and a graduation of 100 grams, the seca 952 fulfills all basic medical standards. Added to that are the typical seca advantages such as fold-up armrests and footrests, secure wheel brakes, an ergonomically shaped seat and smooth-running rollers, which, together with the ergonomically shaped handle, make it very easy to push the chair scale.

## Simple operation.

Thanks to the fold-up armrests and footrests, elderly and heavy patients can take a seat on the scale easily and comfortably. It is just as easy for personnel to use the scale because the batteries eliminate the need for connection to an electrical outlet and keep on working through up to 4000 weighings. No time is lost to a learning curve. A user can operate the scale intuitively. Display elements are conveniently located in the handle and can be set to show kilograms, pounds or stones.

## TARE and HOLD.

Of course the seca 952 is also equipped with the important TARE function which allows every newly added weight to be ascertained separately and the HOLD function which lets medical personnel take care of the patient first and read out the results later. All these advantages plus well-known workmanship and durability make the scale a genuine seca product. It fulfills all the expectations people have of the world leader in medical weighing and measuring.



*The new ergonomic handle makes scale easier to push.*



*The practically positioned display is easy to use.*



*To make sitting down and standing up easier, the armrests and footrests fold up out of the way.*

**seca 952**

## Technical Data

- Capacity: 200 kg / 440 lbs / 32 sts
- Graduation: 100 g / 0.2 lbs
- Dimensions (WxHxD): 565 x 920 x 810 mm
- Device weight: approx. 24.5 kg
- Power supply: batteries
- Functions: TARE, HOLD, kg/lbs/sts switch-over

**ECOMED**

Serving The HealthCare Industry

**T 1800 67 5432**  
**F 1800 67 5433**  
**E enquiries@ecomед.com.au**  
**W ecomед.com.au**